



Image from jax.org

# Mouse Strain Datasheet

## Strain Name

Pkd2<sup>fl/fl</sup>

**MGI Gene ID**

MGI:4843126

**Full Allele Name**

Pkd2<sup><tm1.1Tjw></sup>

**Type of Allele**

conditional allele

**Human Gene (HGNC)**

PKD2

**Genetic Background**

C57BL/6J

**Commercial Source**

Jax mice

**Stock Number**


17292

**Link**

<https://www.jax.org/strain/017292>

**Genotyping Protocol**

<https://www.jax.org/Protocol?stockNumber=017292&protocolID=24847>



## Strain Details

A loxP site was inserted upstream of exon 11 and an FRT-flanked neo cassette with a 5' loxP site was inserted downstream of exon 13. Flp-mediated recombination removed the neo cassette leaving exons 11 through 13 floxed. These Pkd2 conditional mutant mice possess loxP sites flanking exons 11-13 of the polycystic kidney disease 2 (Pkd2) gene. Expression of Cre recombinase results in removal of the intervening exons generating a null allele. This strain may be useful for studying renal development in autosomal dominant polycystic kidney disease or the function of pkd2 in various tissues and throughout adult life.

## Validation or publication

<https://pubmed.ncbi.nlm.nih.gov/20862291/>

**Contact Name**

Patricia Outeda

**Email**

pgarcia@som.umaryland.edu

## PCR Protocol for Genotyping:

Pkd2fl/fl

### A. Digestion of mouse tail or ear notch, and **embryo tail (in red)**:

1. Add 100µL of Tissue Digestion Buffer and 2µL of Proteinase K per tail (~1-2mm length). For embryos tail add **50µL of Tissue Digestion Buffer and 1µL of Proteinase K**. Make sure tail is immersed in the buffer.
2. In a thermocycler incubate at 55°C for 1 h followed by 95°C for 8 min to inactivate the enzyme and hold at 10°C. For embryos incubate at 55°C for 30 min followed by 95°C for 8 min and hold at 10°C.
3. Vortex and store at 4°C (-20°C for long storage) or use immediately to set up the PCR.

### B. PCR Genotyping Protocol

Primers			
flox11-13(A)-F	5'-	CCT TTC CTC TGT GTT CTG GGG AG	-3'
"flox11-13(B)-R	5'-	GTT TGA TGC TTA GCA GAT GAT GGC	-3'
	5'-		-3'
	5'-		-3'

PCR Reaction		PCR Conditions		
BioMix (Bioline)	10.0 µL		Heated Lid	105°C
Primers (@10 µM each)	0.8 µL		Initial Denaturation	94°C 5 min
			Number of Cycles	x35
ddH <sub>2</sub> O	7.2 µL		94°C	20 sec
			56°C	35 sec
			72°C	35 sec
DNA template	2.0 µL		Final Extension	72°C 10 min
Total Volume	20.0 µL		Final Hold	10°C

PCR Product Size (bp)	
Wild type band	232 bp
Pkd2fl band	318 bp

### C. Reagents

Reagent	Cat #	Final Concentration	Working Concentration
<b>Tissue Digestion Buffer for ear notch or tail</b>			
Tris pH8.5		50mM	
EDTA		1mM	
Tween20		0.5%	
<b>Proteinase K (Invitrogen)</b>	25530-015	20mg/mL	
<b>BioMix (Bioline)</b>	BIO-25012		



Image from jax.org

# Mouse Strain Datasheet

## Strain Name

Pkhd1/Cre

MGI Gene ID

n/a

Full Allele Name

Tg(Pkhd1-cre)1lgr

Type of Allele

Cre

Human Gene (HGNC)

n/a

Genetic Background

C57BL/6J

Commercial Source

UTSW O'Brien Kidney Center

Stock Number


Pkhd1Cre

Link

[https://utsw.corefacilities.org/service\\_center/show\\_external/3550?name=utsw-o-brien-kidney-center-animal-models-core](https://utsw.corefacilities.org/service_center/show_external/3550?name=utsw-o-brien-kidney-center-animal-models-core)

Genotyping Protocol

<https://pubmed.ncbi.nlm.nih.gov/18263604/>



## Strain Details

Pkhd1/Cre mice mediate Cre/loxP recombination in late embryonic/postnatal renal collecting ducts, the ureter, and intrahepatic bile ducts.

## Validation or publication

<https://pubmed.ncbi.nlm.nih.gov/18263604/>

Contact Name

Stephen Parnell, Patricia Outeda

Email

sparnell@kumc.edu, pgarcia@som.umaryland.edu

## PCR Protocol for Genotyping: Pkhd1/Cre

### A. Digestion of mouse tail or ear notch, and **embryo tail (in red)**:

1. Add 100µL of Tissue Digestion Buffer and 2µL of Proteinase K per tail (~1-2mm length). For embryos tail add **50µL of Tissue Digestion Buffer and 1µL of Proteinase K**. Make sure tail is immersed in the buffer.
2. In a thermocycler incubate at 55°C for 1 h followed by 95°C for 8 min to inactivate the enzyme and hold at 10°C. For embryos incubate at 55°C for 30 min followed by 95°C for 8 min and hold at 10°C.
3. Vortex and store at 4°C (-20°C for long storage) or use immediately to set up the PCR.

### B. PCR Genotyping Protocol

Primers			
Pkhd1 Forward	5'-	TCTGTCTCAACATAACTCATTG	-3'
Cre Reverse (this is a generic Cre primer)	5'-	TCGACCAGTTTAGTTACCC	-3'
	5'-		-3'
	5'-		-3'

PCR Reaction		PCR Conditions		
BioMix (Bioline)	10.0 µL		Heated Lid	105°C
Primers (@10 µM each)	0.8 µL		Initial Denaturation	94°C 5 min
			Number of Cycles	x35
ddH <sub>2</sub> O	7.2 µL		94°C	20 sec
			60°C	35 sec
			72°C	35 sec
DNA template	2.0 µL		Final Extension	72°C 10 min
Total Volume	20.0 µL		Final Hold	10°C

PCR Product Size (bp)	
Wild type band	NA
Cre (specific for Pkhd1-Cre)	~800 bp

### C. Reagents

Reagent	Cat #	Final Concentration	Working Concentration
<b>Tissue Digestion Buffer for ear notch or tail</b>			
Tris pH8.5		50mM	
EDTA		1mM	
Tween20		0.5%	
<b>Proteinase K (Invitrogen)</b>	25530-015	20mg/mL	
<b>BioMix (Bioline)</b>	BIO-25012		